

Repeat at-home sampling of gamified behavioural tasks and wearable dry-EEG can match the quality of lab-based systems

Barbey FM^{1,2}, Buick AR³, Danyeli L^{4,5}, Dyer JF^{1,2}, Md Islam¹, Krylova M^{4,5}, Murphy B¹, Nolan H¹, Rueda-Delgado LM^{1,2}, Walter M^{4,5}, and Whelan R²

¹Cumulus Neuroscience Ltd., Dublin, Ireland

²Trinity College Dublin, Dublin, Ireland

³Cumulus Neuroscience Ltd., Belfast, United Kingdom

⁴Jena University Hospital, Jena, Germany

⁵Leibniz Institute for Neurobiology, Magdeburg, Germany

Contact: florentine@cumulusneuro.com

Motivation

- Recent advances in sensor development, signal processing, and cloud-based technologies have led to the creation of wireless, dry EEG recording systems (1)
- Gamification of laboratory cognitive paradigms – adding game-like features (points, graphics, levels, storyline, etc.) – reduce attrition and increase participant engagement over extended periods (2,3)
- These technologies combined allows the repeated collection of task-driven EEG data by research participants themselves over days, weeks, or months (4).
- The quality of remote self-administered gamified EEG needs to be investigated

Method

- 30 healthy male adults (mean age=25.6 years) completed 4 weeks of self-administered dry EEG at-home recordings
- Gamified visual Oddball and Flanker tasks
- Usability was evaluated via participant adherence, percentage of sessions successfully completed, and quantitative feedback using the System Usability Scale.
- EEG signals were pre-processed with frequency domain filtering and automatic artefact removal methods. Low quality epochs were removed using a conservative correlational approach, discarding ~50% of trials.

- Signal variability was quantified by computing the Standardized Measurement Error (SME) (5) over repeated sessions and multiple trials. Event-related potentials (ERPs) from the Oddball task were extracted from stimulus locked Target trials at Pz. SMEs were calculated across the 400–500 ms time window. ERPs from the Flanker task were extracted from response locked Correct and Error trials at FCz. SMEs were calculated across the 0–100ms time window. The shaded areas correspond to the 95% confidence intervals.

- Wet EEG reference SME values were extracted from the open-source ERP CORE EEG database for comparison (6).



ERPs analyses

Top Row: Event-related potentials. **Middle row:** Standardized Measurement Errors per number of trials. The grey continuous lines correspond to the number of participants remaining in the SME calculation after trial rejection. **Bottom Row:** Mean numbers of trials available after preprocessing. The black dotted lines correspond to the mean numbers of trials extracted from the wet EEG study.




Results

- 684 EEG sessions were collected. On average, participants completed 22.87 sessions of 23 sessions requested (std=1.92).
- Participants successfully completed 96% (std=5%) of the sessions attempted and reported a mean usability score of 88.3 (std=9.90).
- In the Oddball, 33 Target trials (red line in the figure), were required to reach the same SME values obtained with 30.5 trials using a wet EEG system.
- For the Flanker task, 467 Correct trials and 56 Error trials were required to reach the same SME values obtained with 337 Correct trials and 40 Error trials using a wet EEG system.
- Allowing for differing epoch rejection rates, aggregation of ERPs across 2 to 4 sessions (depending on task) resulted in grand average signal quality comparable to single-session data collected using wet EEG.

Cumulus Headset



Flanker and Oddball Gamified Tasks



Conclusions

Wireless EEG technology is a suitable tool for cognitive neuroscience investigations, and has the potential to provide objective, frequent and patient-centered tracking of biomarkers of functional neurophysiology. This approach has potential to facilitate large scale longitudinal studies of brain disorders that manifest on different time scales.

Bibliography

- Lau-Zhu A et al, Dev Cogn Neurosci (2019)
- Germine L et al, Neuropsychopharmacology (2020)
- Lumsden J et al, JMIR Serious Games (2016)
- McWilliams EC et al, Front Psychiatry (2021)
- Luck SJ et al, Psychophysiology (2021)
- Kappenman ES et al. NeuroImage (2021)

